

TMT™ TOURNIQUET**\$29.95**

Tighten With Confidence

**SKU:** 31-115**NSN:** 6515-01-656-6191 | **Info::** Made in USA |**PRODUCT DESCRIPTION****CoTCCC Recommended**

The Tactical Mechanical Tourniquet (TMT™) is a hemorrhage control device specifically designed for massive hemorrhage control of an extremity. Using lessons learned from the battlefield, the TMT is designed to ensure ease of application and definitively achieve occlusion. Taught as a self-aid and buddy-aid task, the TMT requires minimal training to instantly treat life-threatening hemorrhage of an extremity. Additionally, the TMT eliminates shortcomings found in other tourniquets.

PRODUCT INFO

- **FAST:** True One-Handed Application in 30 seconds or less
- **DURABLE:** Co-poly resin acetyl with IR/UV protective additive
- **BETTER:** eliminates the current tourniquet “percentage of failure” rate

Product Attributes

- True one-handed application for both upper and lower extremities in 30 seconds or less
- Wider band for occlusion at lower pressure and increased comfort
- Dual locking mechanism to prevent slippage and loss of pressure
- Audible signature “click” identifies torsion bar is secure
- Advanced materials with IR and UV protection
- Integrated UV protection for increased durability
- Accommodates varying limb circumferences:
 - 5-3/4” Min

- 38" Max

Clinical Benefits

- Allows for faster application in one handed and two handed application compared to other tourniquets
- Increased comfort from the wider band allows for better occlusion at lower pressure without further harm to the patient
- Ease of application and intuitive use
- Cost effective-
 - Reduced budget requirement for fielding of force
 - Increased capability and efficacy over current solution

SPECIFICATIONS

- NSN: 6515-01-656-6191
- PN: 31-115
- CoTCCC Recommended
- FDA Listed / CE Marked
- Not made with natural rubber latex / Single patient use
- Shelf life:
 - Indefinite stored shelf life (controlled environment)
 - Up to 10 year shelf life (operational storage)
- AbilityOne manufactured / contract
- Dimensions-
 - Packaged: 4.5"L x 2.5"W x 1.85"D
 - Open: 38.7"L x 2.25"W
 - Weight: 2.9 oz.
 - Accommodates limbs:
 - Min: 5 3/4"
 - Max: 38"
- Co-poly resin acetyl with IR/UV protective additive
- US Government owned design
- Berry Compliant
- Made in USA

TRAINING

[NAMRU-SA Report 2015-52](#)

[NAMRU-SA Report 2014-23](#)

[TMT IFU](#)

[Phase IIIA Study Report](#)

[Phase IIIB Study Report \(Redacted\)](#)

[Tourniquet Video](#)