

TRI-FIT SILICONE EARPLUGS

\$10.64

Preserve Your Hearing



SKU: N/A

NSN: Small - NSN 6515-00-442-4821,
Medium - 6515-00-442-4818, Large - NSN
6515-00-467-0092 | **Info::** Pk of 12 Pair |

PRODUCT DESCRIPTION

Tri-Fit Silicone Earplugs are made of soft, pre-molded, hypo-allergenic silicone offering unsurpassed comfort and protection. Sizes are conveniently color-coded (small- green, medium - orange, blue - large) for easy identification.

PRODUCT INFO

- NRR=25db
- Three ultra thin flanges provide exceptional fit while minimizing pressure to the inner ear
- May be reused indefinitely
- Washes in warm, soapy water without compromising flexibility, comfort or protective capabilities
- Color coded by size
- Inserter aides in proper insertion
- Tested in an independent accredited NVLAP facility

SPECIFICATIONS

- **PN: 40-201 (small, green)**
- NSN: 6515-00-442-4821
- **PN: 40-202 (medium, orange)**
- NSN: 6515-00-442-4818
- **PN: 40-203 (large, blue)**
- NSN: 6515-00-467-0092
- Unit weight (PK /12 pair): 0.8 oz

- Made in USA