



TRUECLOT® HEMORRHAGE TRAINING KIT

\$625.00

Realistic blood for realistic training

SKU: 30-561

PRODUCT DESCRIPTION

TrueClot® Hemorrhage Training Kit (Large Caliber Gunshot Wound) provides realistic, low-cost and repetitive training for wound packing with the use of advanced hemostatic gauzes. The task trainers have very realistic outer anatomy and wound beds with hidden bleed points. They are hand made from sturdy silicone and contain no moving parts. The Task Trainer also has a bony landmark that allows for deep packing instruction for large caliber gunshot wounds. TrueClot Blood Simulant is non-biological, closely mimics human blood, and has a 12 month shelf life. It is non-staining and easily washable from skin and clothing. The entire kit contains: TrueClot Blood Simulant Concentrate Packs (4), TrueClot Wound Packing Task Trainer, TrueClot Simulated Hemostatic Gauze (9), mixing bottle in protective hard case.

PRODUCT INFO

- PN: 30-561
- Gov't Open Market Price: \$625.00
- Cost-effective – contents can train approx. 20 people
- Realistic hemorrhage control for TCCC, TECC, TEMS, PHTLS, Active Shooter Training
- TrueClot® Blood Simulant
 - 12 month shelf life
 - Non-biological, safe to use
 - Closely mimics human blood
 - Washes from skin and clothing with soap & water

SPECIFICATIONS

- Task Trainer Dimensions: 7.25" L x 6" W x 2.5" H
- Storage Case Dimensions: 13.5" x 10" x 5.5"
- Liquid Capacity: 1 Liter