

**UNIVERSAL SPLINT SET****\$18.14**

---

**SKU:** 55-120**NSN:** 6515-01-688-6388 |**PRODUCT DESCRIPTION**

The Universal Splint Set is designed for emergency management of bone fractures. Proven for treating long bone fractures, it can be used for the humerus, lower arm, wrist, foot and lower leg. Its exclusive flame retardant, textured foam provides an additional safety feature for the warfighter. Once the splint is properly placed, it can be fastened using the enclosed, muslin securing ties. The waterproof, textured foam has both antimicrobial and antifungal properties, and the pebble skin finish reduces skin tear on the foam while allowing for less slippage when moving the patient. The entire perimeter of the splint is fusion bonded, eliminating the need for adhesives which can breakdown over time.

**PRODUCT INFO**

- Flame Retardant
- Includes 4 Muslin Securing Ties
- Antibacterial/Antifungal Properties
- Malleable, Aluminum Alloy Core
- Lightweight/Low Cube
- Complies with MIL STD 810F
- One Size Fits All
- Radiolucent
- Waterproof
- Reusable

## **SPECIFICATIONS**

- Non Shelf Life Product
- Specifications:
  - Dimensions:
    - Packaged: 4" x 9" x 1"
    - Deployed Splint: 4" x 36"
    - Securing Tie: 2" x 36" (EA)
  - Packaged Weight: 5 oz
- Coyote Brown Color
- Latex Free
- Non sterile
- Berry Compliant
- Made in USA