



WARRIOR BLOOD & FLUID WARMER SET (STANDARD SET)

\$4,449.00

No More Cold Blood

SKU: 37-600

PRODUCT DESCRIPTION

Recognized for its superiority over prevailing warmers in the recently released CoTCCC ARC Guidelines, the Warrior Blood & Fluid Warmer is a high performance modular system that uniquely meets blood and IV fluid warming needs across the entire continuum of emergency care, simplifying patient handoff and reducing costs.

With its use of medical-grade stainless steel for warming, the Warrior is the only portable solution capable of warming fluids from practically any fluid input temperature to 38°C in real time, even at high flow rates. The device's unmatched capacity provides three to five times the amount of warmed fluids per single battery than any alternative solution.

The complete set contains the Warrior Standard Base Unit, Compact Disposable Unit (2), Rechargeable Enhanced Battery, Battery Charger and Adapter in a Soft Carry Case.

PRODUCT INFO

- Warms near frozen fluid to 38°C in under 11 seconds
- Disposable Unit provides unlimited volume to single patient
- Up to 3 liters at 4° C fluid input temperature / up to 5 liters at 20°C fluid input temperature using single battery
- Functions on gravity feed
- 400+ charging cycles per battery
- Two disposable unit (DU) options (Standard or Compact)
- High flow rates 200/290 ml/mn (battery/AC) for full warming range (4-38°C/39-100.4°F)
- Does not require proprietary IV line

- Reduced lifecycle cost when compared to Thermal Angel and Belmont Buddy1
- Simple to operate
- Portable: Can fit small transport platforms
- Highly efficient technology: 3-5 liters of warmed fluids with a single battery
- Mountable to pole, rail or stretcher
- Communicative, built-in display
- No need for periodic calibration
- Affordable consumables
- Practically Zero maintenance: 5 years between service cycles
- Patent-Protected Technology
- Field Proven Technology
- Same consumable can be used across the entire continuum of care

SPECIFICATIONS

- PN: 37-600 (Soft Case Set- Standard Configuration)
- FDA 510(k) Cleared
- Base Unit:
 - 5 Year service life
- Disposable Unit:
 - Single patient use
 - 3 Year shelf life
- Specifications:
 - Set Dimensions: 16"L x 12"H x 6.5"W
 - Set Weight: 8 lbs (approx.)