



WARRIOR EXTREME BLOOD & FLUID WARMER SET

\$5,850.00

EXTREME Military enhancements

SKU: 37-604

PRODUCT DESCRIPTION

Recognized for its superiority over prevailing warmers in the recently released CoTCCC ARC Guidelines, the Warrior EXTREME is a high performance warming system that uniquely meets blood and IV fluid warming needs across the continuum of emergency care, simplifying patient handoff and reducing costs.

Using medical-grade stainless steel for warming, the Warrior is the only portable solution capable of warming fluids from practically any fluid input temperature to 38°C in real time, even at high flow rates. The device's unmatched capacity provides three to five times the amount of warmed fluids per single battery than alternative solutions.

The Warrior EXTREME has the same form factor as the current, standard Warrior version, but with several enhancements designed for the military market, to include: hardened circuitry for air worthiness, Trade and Berry compliant. The Warrior EXTREME is certified for the latest, most advanced EMS standard, IEC 60601-1-12 (and IEC 60601-1-2, 4th edition). It is 3X more immune to radiated electromagnetic fields, 2X more immune to electrostatic discharge, 10X more immune to magnetic field, and has enhanced ingress protection (IP) rating (IP33 or higher) as well as better resiliency to drop, impact and vibration.

The complete set contains the Warrior EXTREME Base Unit, Compact Disposable Unit (2), Rechargeable Battery (2), Battery Charger, Adapter and Instructions for Use in a Berry Compliant, Soft Carry Case.

PRODUCT INFO

- Warms near frozen fluid to 38°C in under 11 seconds
- High flow rates up to 200 ml/mn
- Disposable Unit provides unlimited volume for a single patient
- Up to 3 liters at 4° C fluid input temperature per battery
- Functions on gravity feed
- 400+ charging cycles per battery
- Does not require proprietary IV line
- Improved performance and reduced lifecycle cost when compared to current solutions

- Simple to operate
- Portable: Can fit small transport platforms
- Highly efficient technology: 3-5 liters of warmed fluids with a single battery
- Mountable to pole, rail or stretcher
- Communicative, built-in display
- No need for periodic calibration
- Practically Zero maintenance: 5 years between service cycles
- Patent-Protected Technology
- Field Proven Technology
- Same consumable can be used across the entire continuum of care

SPECIFICATIONS

- PN: 37-604 (Soft Case Set- EXTREME Configuration)
- FDA 510(k) Cleared
- Ruggedized, hardened circuitry for air worthiness
- Berry Compliant
- 100% TAA Compliant
- Base Unit:
 - 5 Year service life
- Disposable Unit:
 - Single patient use
 - 3 Year shelf life
- Specifications:
 - Set Dimensions: 16"L x 12"H x 6.5"W
 - Set Weight: 8 lbs (approx.)