

WATER-JEL BURN BLANKETS

From: \$380.00

SKU: N/A



PRODUCT DESCRIPTION

Water-Jel Fire Blankets are a one-step system for emergency first aid on a burn victim. The blanket puts out the fire, stops burn progression and cools the skin, relieves pain and protects the burn area from contamination.

Water-Jel Fire Blankets consist of a scientifically formulated cooling gel and a special carrier. The 100% worsted wool blanket with intercellular weave can absorb up to 13 times its own weight. The blanket is soaked in a water-based / water-soluble, bacteriostatic and biodegradable gel. Can be placed directly over clothing and won't stick to burn injuries. Shelf life: 5 years from date of manufacture.

SPECIFICATIONS

3'X2.5" Burn Wrap:

Packaged in ruggedized coyote brown twist top canister.

Case Pack: 4 canisters individually boxed

Case Weight: 26.5 lbs.

Case Size: 15.5"Hx 13.5"W x 6.75"D

Blanket Dimensions/Unit: 3' x 2.5' (91 x 76 cm)

5'X6' Burn Blanket:

Packaged inside a nylon coyote brown pouch.

Case Pack: 4 pouches individually boxed

Case Weight: 56 lbs.

Case Size: 19.5"H x 11.75"W x 11"D

Blanket Dimensions / Unit: 6' x 5' (183 cm x 152 cm)

5'X6' Burn Blanket:

Packaged in ruggedized coyote brown twist top canister.

Case Pack: 4 canisters individually boxed

Case Weight: 82 lbs.

Case Size: 18.5"H x 15.75"W x 18"D

Blanket Dimensions / Unit: 6' x 5' (183 cm x 152 cm)