

**XTRACT® 2 LITTER****\$369.00**

For Extreme TACEVAC

SKU: 52-140**PRODUCT DESCRIPTION****The Ideal Solution for Rapid Casualty Extraction**

In extreme environments, the Xtract® 2 Litter is the tactical rescue solution for evacuating casualties safely and efficiently. It is the ideal product for responders who require a rugged, flexible, low cube, litter. Casualties can be quickly dragged or carried from point of injury to safety.

Developed with input from the battlefield, the Xtract2 cocoons the casualty within the fabric, protecting existing treatments such as tourniquets and IV lines. The waterproof, heat-retaining fabric is tear resistant and includes handles and securing straps strategically-placed to support the casualty, maintain alignment and reduce secondary injury.

PRODUCT INFO

- Secures & Protects the Casualty
- Light Weight and Functionally Robust
- Reduces Secondary Injury

- Tear-resistant, double-layered fabric
- Rated to carry 600lbs
- Fire retardant inner liner
- Simple, intuitive use for all skill levels
- Accommodates a sleeping mat between fabrics to increase insulation, rigidity and comfort
- Cocoons patient within fabric, protecting existing treatments such as tourniquets and IV lines

- Minimizes patient movement
- Reduces risk of secondary injury with strategic placement of handles and securing straps for optimal anatomical support
- Xtract2 fabric pattern helps prevent conventional heat loss, maximizing additional hypothermia equipment
- Accommodates different casualty lengths
- Meets the required 5th to 95th percentile population range

SPECIFICATIONS

- Unit Weight: 3 lbs
- Packed: 11.5" x 5.5"
- Length (unadjusted): 7' 6"
- Length (adjusted): 4' 8"
- Maximum load: 660 lbs

TRAINING

[Xtract®2 adopted by British Military](#)

[Xtract Training Video](#)