

Testing Parameters:

Stainless steel  
15 minutes  
.5 second delay  
OZ  
20 second test  
12 in./ second

Test	Average (oz)
T1	107.6
T2	132.42
T3	127.01
T4	122.18
<b>Average force</b>	<b>122.3025</b>

Test	Average (oz)
T1	37.59
T2	37.22
T3	40.5
T4	37.15
<b>Average Force</b>	<b>38.115</b>

Testing Parameters:

Stainless steel  
1 minute  
.5 second delay  
OZ  
20 second test  
12 in./ second

Test	Average (oz)
T1	122.12
T2	99.79
T3	108.88
T4	121.18
<b>Average Force</b>	<b>112.9925</b>

Test	Average (oz)
T1	32.85
T2	41.88
T3	40.01
T4	38.97
<b>Average Force</b>	<b>38.4275</b>

Testing Parameters:

Skin  
15 minutes  
.5 second delay  
OZ  
20 second test  
12 in./ second

Test	Average (oz)
T1	10.99
T2	14.92
T3	15.51
T4	19.04
<b>Average Force</b>	<b>15.115</b>

Test	Average (oz)
T1	6.26
T2	8.47
T3	6.66
T4	11.77
<b>Average Force</b>	<b>8.29</b>

Testing Parameters:

Skin  
1 minute  
.5 second delay  
oz  
20 second test  
12 in./ second

Sentinel Chest Seal	
Test	Average (oz)
T1	22.72
T2	24.14
T3	27.43
T4	30.42
<b>Average Force</b>	<b>26.1775</b>

Asherman Chest Seal	
Test	Average (oz)
T1	6.75
T2	9.62
T3	7.13
T4	9.61
<b>Average Force</b>	<b>8.2775</b>