

Testing Parameters:

Stainless steel
15 minutes
.5 second delay
OZ
20 second test
12 in./ second

Test	Average (oz)
T1	107.6
T2	132.42
T3	127.01
T4	122.18
Average force	122.3025

Test	Average (oz)
T1	37.59
T2	37.22
T3	40.5
T4	37.15
Average Force	38.115

Testing Parameters:

Stainless steel
1 minute
.5 second delay
OZ
20 second test
12 in./ second

Test	Average (oz)
T1	122.12
T2	99.79
T3	108.88
T4	121.18
Average Force	112.9925

Test	Average (oz)
T1	32.85
T2	41.88
T3	40.01
T4	38.97
Average Force	38.4275

Testing Parameters:

Skin
15 minutes
.5 second delay
OZ
20 second test
12 in./ second

Test	Average (oz)
T1	10.99
T2	14.92
T3	15.51
T4	19.04
Average Force	15.115

Test	Average (oz)
T1	6.26
T2	8.47
T3	6.66
T4	11.77
Average Force	8.29

Testing Parameters:

Skin
1 minute
.5 second delay
oz
20 second test
12 in./ second

Sentinel Chest Seal	
Test	Average (oz)
T1	22.72
T2	24.14
T3	27.43
T4	30.42
Average Force	26.1775

Asherman Chest Seal	
Test	Average (oz)
T1	6.75
T2	9.62
T3	7.13
T4	9.61
Average Force	8.2775