



HEAD & HYPOTHERMIA





EYE SHIELD

1 ea - 2.5" x 3.25", aluminum eye shield w/garter




Step 1:  ID wound. Perform a rapid, field visual acuity test.



Step 2: Cover injured eye with eye shield. Secure as shown. Reassess. 



 **Combat Medical**
5555 Harrisburg Industrial Park Drive
Harrisburg, NC 28075 USA
1-855-428-6074
www.combatmedicalsystems.com

REF 06-200 CMS-PUB-EYE-04-004



Made in USA

WARNING: For use by medical providers. Place no objects between eye and shield.

Pat. Info: combatmedicalsystems.com/ip
Prod. Info: combatmedicalsystems.com/e-labels
©2018 Combat Medical Systems, LLC