



HEAD & HYPOTHERMIA





EYE SHIELD

1ea - 2.375" x 3", polycarbonate eye shield w/garter




Step 1:  ID wound. Perform a rapid, field visual acuity test.



Step 2: Cover injured eye with eye shield. Secure as shown. Reassess. 



 Combat Medical
5555 Harrisburg Industrial Park Drive
Harrisburg, NC 28075 USA
1-855-428-6074
www.combatmedicalsistemas.com



Made in USA

REF 06-205 CMS-PUB-PCEYE-04-003

WARNING: For use by medical providers. Place no objects between eye and shield.

Pat. Info: combatmedicalsistemas.com/ip
© 2017 Combat Medical Systems, LLC